

1 Securely attach one end of the tubing to the device.



2 Securely attach the other end of the tubing to the flow meter.




3 Set the desired flow rate (starting at 5 L/min), slowly increasing to reach the desired expiratory airway pressure. Do not exceed a flow of 15 L/min.



4 As the patient exhales through the device, verify that the manometer gauge is at the desired pressure range.



5 Have the patient inhale and exhale for the prescribed amount of time.



6 Repeat steps 1 through 4 as many times as instructed by your physician or health care professional.

These are abbreviated instructions for use. Please read the complete instructions that are packaged with the VersaPAP™ device.