





reliable oximetry


Philips Nasal Alar SpO₂
Sensor by Curaplex®





the ala: highly perfused and easily accessible

The alar sensor measures blood oxygen levels at the ala – where the nose meets the cheek and facial arteries are fed by the internal and external carotid arteries – to produce a strong, consistent signal. In contrast, when a patient is in shock or when vasopressors limit the blood supply to the extremities, finger sensors can become unreliable. The ala is also easily accessible during surgeries in which patients are supine, and the nose moves very little in comparison to the digits.



reliable results for confident assessment

Measuring arterial oxygen saturation in poorly perfused patients often results in wasted time and discarded sensors as you search for a strong signal. Yet for these patients, accurate and responsive measurement is critical. The **Philips Nasal Alar SpO₂ Sensor by Curaplex®** measures oxygenation at the ala, which is an optimal measurement site for these difficult-to-measure patients. So you can be confident in your readings, even during critical states of low perfusion and blood flow centralization.¹

experience the benefits of alar sensing

- Patient's hands are sensor-free
- Sensor design aids in correct orientation application
- Easy accessibility during procedures
- One sensor for extended hospital stays
- Adhesive free
- Lower occurrence of pressure injury than forehead sensors¹
- Wipeable

comfortable and durable

Trying – and discarding – multiple sensors when seeking a steady signal can be frustrating and costly. The alar sensor is easy to place and delivers a dependable signal right from the start. With a durable clip and molded, medical-grade silicone, the alar sensor is designed for comfort. There is no adhesive, limited pressure, little heat and no annoying headband. You simply move the sensor to the opposite ala as instructed, and you can use the same sensor for up to seven days, which easily outlasts most hospital stays.



critical care

The perfusion of hemodynamically unstable patients results in caregivers having to spend excessive time trying to find a stable SpO₂ signal. The search for a quality signal may include multiple site rotations and frequent trips to the stockroom. Even in hemodynamically unstable patients, such as patients on vasopressors or in shock states, the ala of the nose is more likely to be well perfused.



surgery

While SpO₂ monitoring has value throughout the hospital, it is a particularly critical measure for anesthesiologists assessing patients during surgery. The alar sensor is responsive to changes in oxygen saturation. After surgery, patients can be moved from the operating room to the PACU, the ICU and beyond, without changing the sensor.



acute care

With mobile patients and heavy demands on caregivers, the alar sensor keeps patient's hands sensor-free so they can more easily perform daily activities. The alar sensor is quick and easy to apply and rotate. Patients can perform frequent hand hygiene unencumbered by typical SpO₂ sensors, reducing their risk for hospital-acquired infection.²

1. Shallom M, Prentice D, Sona C, Mazuski J. Comparison of Nasal and Forehead Oximetry Accuracy and Pressure Injury in Critically Ill Patients. Critical Care Medicine. 2016;44:12(Suppl.).
2. Haverstick S, Goodrich C, Freeman R, James S, Kullar R, Ahrens M. Patients' Hand Washing and Reducing Hospital-Acquired Infection. Critical Care Nurse. 2017;37(3):e1-e8.

the nasal alar can be a hospital-wide solution for oximetry

low perfusion solution components

plug-in module for monitor*

description	qty	item #
Philips FAST SpO ₂ Module	1/ea	M3-M1020B

*Monitor sold separately, contact us for more information



alar sensor

description	technology	single patient use	sterile	application site	patient application	latex free	cable length	shelf life	qty	item #
Nasal Alar SpO ₂ Sensor	FAST	Yes, up to 7 days	No	Nasal ala(r)	Pediatric and Adult	Yes	3.28 ft	5 Years	20/ca	73-989803205381
	Multi-compatible								24/ca	301-11214

cable connection

description	type	cable length	qty	item #
Nasal Alar Oximetry Cable	8-pin to 9-pin adapter	9.8 ft	1/ea	73-M1943AL



Contact us to walk through your recommended oximetry configuration.



Ask us about these items and our Rewards Program!

tri-anim.com/thc/capital-rewards