



AccuPAP™

The AccuPAP™ Positive Airway Pressure device is an innovative treatment option used to provide positive airway pressure during lung volume expansion therapy. Indicated for the prevention & treatment of atelectasis, AccuPAP can aid in reducing the air trapping in asthma and/or COPD patients as well as optimize delivery of bronchodilator medications.

- Adjustable pressure valve eliminates the need for manometer
- Treatment can be delivered with mouthpiece or mask attachment
- Design allows for convenient delivery of aerosol medication during treatment
- Pressure port with cap and ambient air inlet
- One step connection to gauge tubing and flow meter
- Available with or without manometer



ORDERING INFORMATION

301-6000	AccuPAP Positive Airway Pressure Device, Mouthpiece, 7' Tubing
301-6001	AccuPAP Positive Airway Pressure Device, Mouthpiece, 7' Tubing, with Manometer
301-5558	Disposable Pressure Manometer