

Eye Protectors



Curaplex® eye protectors are rigid, clear, plastic protective covers used to shield the patients' eyes during surgical procedures. The eye protectors feature a non-allergenic, self-adhesive foam cushion to enable quick and accurate application on the patient. The thick foam padding increases patient comfort even when the patient is in prone position. The one-piece backing has "courtesy tabs" to assist with removal. Curaplex® eye protectors are latex free, sterile and single patient use. Made in the USA.

Ordering Information

Adult Eye Protector, 25 ea/bx
301-AEG100

For more information, please contact your Account Manager or call Customer Service at 800-TRI-ANIM (874-2646).



Tri-anim®

800.874.2646 • www.tri-anim.com