

# Anesthesia Breathing Circuits

**Secure Connections** are welded together to prevent leakage and maintain connection.

**High Quality Tubing** made of stronger material to minimize risk of leakage with excellent drapability.

**Secure Hand Grip** enables easy removal from the machine.

**Tethered Gas Sampling Ports** are uniquely attached to prevent misplacement or loss.

**Soft, Tactile Breathing Bag** provides excellent feel with a secure connection to the machine.

**Soft Cushion Anesthesia Masks** provide superior seal and patient comfort.



Available in parallel and single limb as well as corrugated and expandable.

Available in both adult and pediatric sizes.

Available in all lengths.

Full range of associated anesthesia filters, connectors, and tubing.

Infinite customization.