

3-in-1 Tubing Adapter

- Quiet; Minimal patient stimulation during use
- Connects Large Bore (22mm), Small Bore (15mm) & 1/4" Cannula to any DISS fitting
- 2-in-1 Option for High Flow applications available (22mm & 15mm)
- Aluminum option coming soon
- Works great with the MaxFLO₂



About the 3-in-1 Adapter

The 3 in 1 adapter is designed specifically for use with multiple sizes of commonly used medical tubing. It gives the user the option to connect 22mm Large Bore, 15mm or 1/4" cannula tubing to a device with DISS fittings. A 2-in-1 style is available for High Flow applications (22mm and 15mm connection).

Specifications

Threads:	Oxygen DISS
Material:	Delrin (white) or Aluminum
Outer diameter:	22mm Male
Inner Diameter:	15mm Femal
Cannula Connection:	1/4" Tapered Barb
Weight:	17 grams (delrin) or 33 grams (aluminum)

Part Numbers

3-in-1 Adapter (1 ea)	R219P31	2-in-1 Adapter (1 ea)	R219P32
Aluminum 3-in-1 Adapter (1 ea)	Coming Soon	Aluminum 2-in-1 Adapter (1 ea)	Coming Soon