

AeroEclipse® XL

Reusable Breath Actuated Nebulizer

OMBRA™

Compressors

ComfortSeal® Masks Available



- Masks contour gently to the face providing a secure and comfortable fit
- Does not contain natural rubber latex, soft silicone elastomer to safeguard against facial irritation
- Anatomically designed for minimal dead space

Nebulizer and Accessory

Description	Case Size	Part Number
Reusable AEROECLIPSE® XL Breath Actuated Nebulizer (includes nebulizer with mouthpiece, EZ Twist Tubing and Owner's Manual)	10	69010
ComfortSeal® Mask with Elbow Adapter, Small	10	10550394010
ComfortSeal® Mask with Elbow Adapter, Medium	10	10550294010
ComfortSeal® Mask with Elbow Adapter, Large	10	10550494010
EZ Twist Tubing, Replacement	10	66510

Compressor Systems

Description	Case Size	Part Number
Ombra™ 120 V Table Top Compressor (includes two AEROECLIPSE® XL BAN, 5 Replacement Air Filters, Owner's Manual)	1	63506
Ombra™ Compressor Air Filters, Replacement (5/pack)	1	66501

Introducing a complete delivery system for true on-demand therapy



Put your patients in control of their therapy

- Breath actuated delivery only produces aerosol on inhalation allowing for breaks in treatment without drug loss
- Reduces potentially harmful emissions for safer home and hospital environments
- Patient friendly, custom designed compressors

Breath-actuation results in a greater inhaled drug mass and less medication waste.¹

AeroEclipse® XL

Reusable Breath Actuated Nebulizer

OMBRA™

Compressors



monaghan™

Contact Information:

Monaghan Medical Corporation
5 Latour Avenue, Suite 1600
Plattsburgh, NY 12901
Customer Service, 800-833-9653
customerservice@monaghanmed.com
www.monaghanmed.com



¹ Sabato K, Ward P, Hawk W, Gildengorin V and Asselin JM. Randomized Controlled Trial of a Breath-Actuated Nebulizer in Pediatric Asthma Patients in the Emergency Department. Respiratory Care 2011;56(6):761 – 770.

SUPERIOR, ASSURED AEROSOL DELIVERY




AeroEclipse® XL BAN
is the *only* breath actuated
jet nebulizer

Breath Actuated Delivery

- Improves respirable dose for efficient aerosol delivery deep into the lungs
- Minimizes environmental loss by producing aerosol only on inhalation
- Offers efficient treatment times to help increase patient compliance

	Respirable Dose	Respirable Fraction (<4.7 µm)	MMAD†
AEROECLIPSE® XL BAN	447.1 µg	70.8%	3.4 µm
LC Sprint†	403.9 µg	57.9%	4.1 µm
SideStream Plus†	208.3 µg	68.6%	3.4 µm

TMI Aerosol Lab, Data on File. 2.5ml, 2.5 mg salbutamol sulphate (Vt = 600 mL, i:E 1:2, 10 bpm) @ 8 minutes † Mass Median Aerodynamic Diameter



OMBRA™ Compressors
Designed for optimal performance with
AEROECLIPSE® XL Breath Actuated Nebulizer



Ombra™ 120 Volt Table Top Compressor

- Sturdy and easy to use for home, clinic or hospital environments
- Average treatment times of 6–8 minutes
- Integrated handle with nebulizer cup holder
- Includes 2 **AEROECLIPSE® XL BAN**

Maximum Pressure	43.5 psi
Free Air Flow Rate	10.0 L/min
Weight	3.5 lbs
Dimensions	7" × 6.1" × 3.9"

Your Choice for Confidence in Aerosol Delivery

www.monaghanmed.com

