



# The first choice for comfort and fit

**Respironics PerforMax** Full Face Mask for multi-patient use (EE, SE)

The Respironics PerforMax Full Face Mask improves comfort, minimizes skin break-down, and eliminates nasal bridge challenges by sealing around the perimeter of the face, where patients have less pressure sensitivity and smoother facial contours. The PerforMax Mask also simplifies mask selection and speeds setup by fitting a wide range of patients with only two sizes.

#### Key advantages

- Standard adult size fits most patients.
- Perimeter facial contact maximizes comfort and simplifies fitting by sealing around the flat surfaces of the face and avoiding the sensitive nasal bridge.
- The clear, unobstructed view of the PerforMax Mask helps reduce patient claustrophobia.

# PHILIPS

# Respironics PerforMax

The Respironics PerforMax Full Face Mask for multi-patient use has additional advantages to improve patient comfort and noninvasive ventilation (NIV) performance.

- Interchangeable elbows accommodate quick transitions between single- and dual-limb ventilator circuits.
- With a choice of the standard elbow (SE) for dual-limb systems and the entrainment elbow (EE) for single-limb systems, clinicians can use the PerforMax Mask with any noninvasive ventilator.
- Snap clips require minimal effort to attach and remove headgear.
- Radial diffusion leak ports reduce noise and deflect exhaled gas away from clinicians.
- Velcro® grab tabs allow quick, easy removal and adjustment of headgear.
- Digital Auto-Trak compatibility enhances NIV performance.
- Simplify reprocessing with single-use replacement headgear.
- Optional sealing pad minimizes air leaks and localized skin breakdown from nasogastric (NG) tubes.

Description	Quantity	Part Number
PerforMax Full Face Mask EE, small adult	1	1047426
PerforMax Full Face Mask EE, adult	1	1047427
PerforMax Full Face Mask EE/SE, small adult	1	1049037
PerforMax Full Face Mask EE/SE, adult	1	1049038
PerforMax entrainment elbow (EE) replacement kit	1	1045313
PerforMax standard elbow (SE) replacement kit	1	1043210
PerforMax small headgear	9	1045325
PerforMax large headgear	9	1045314
Accessory swivel replacement kit (22mm male-male)	1	7041
Nasogastric tube sealing pad	10	452040



Entrainment elbow (EE)



Standard elbow (SE)



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved. PerforMax and Respironics are trademarks of Respironics, Inc., dba Philips Respironics, and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)  
 fax: +31 40 27 64 887

Hoech JP MCI 4102187 PN 1060288  
 4522 962 43141 \* FEB 2009

Philips Healthcare  
 Global Information Center  
 P.O. Box 1286  
 5602 BG Eindhoven  
 The Netherlands