

SALTER LABS®



Salter-Style® Adult Cannulas

Salter's unique, patented design results in greater patient comfort



Features

- Anatomically curved to fit upper lip
- Curved, tapered nasal prongs
- Over the ear style
- Three-channel safety tubing
- Soft, lightweight, clear material

Benefits

- Eliminates flap or ridge which can cause irritation
- Improved anatomical fit
- Secure positioning
- Uninterrupted flow of oxygen
- Better comfort for long-term use. Does not interfere with patient observation

Wrap around, "Anti-Roll" Feature


Works with the anatomy rather than against it...A Salter Innovation.



Salter-Style® Adult Cannula 1600

Clear, lightweight cannula with unique one piece anatomical design. Eliminates the need for flap, which can cause irritation. The result is maximum patient comfort. Over-the-ear style ensures secure positioning, allowing complete freedom of movement. Provided with three-channel safety tubing or smooth bore tubing in lengths from 4 to 100 feet. Special lengths available.

One piece anatomical design



Three-channel safety tubing

Non-flared tips

Smooth bore or three-channel safety tubing



Adult Conventional Cannula 1056

Conventional cannula with non-flared tips. Made of soft disposable material. Adjustable tubing fits comfortably over the ears to hold the cannula firmly in place. Provided with standard smooth bore tubing in 7' or 16' lengths. Special lengths available. See price list for details.

Ordering Information

Adult Nasal Cannulas*

Salter-Style® Cannulas

| | Units (per case) | Safety Channel | Smooth Bore |
|--|------------------|-------------------|-------------|
| Clear cannula with 7' supply tube | 50 | 1600 ^x | 1650 |
| Clear cannula without supply tube | 50 | 1606 ^x | |
| Clear cannula with 7' supply tube and female thread grip connector | 50 | 1699 | |

Special Length Salter Cannulas*

| | Units | Safety Channel | Smooth Bore |
|--|-------|---------------------|-------------|
| Adult cannula with 1' supply tube | 50 | 1600-1 | |
| Adult cannula with 4' supply tube (ideal for strollers and portable units) | 50 | 1600-4 ^x | 1650-4 |
| Adult cannula with 9' supply tube | 50 | 1600-9 | 1650-9 |
| Adult cannula with 10' supply tube | 50 | 1600-10 | 1650-10 |
| Adult cannula with 12' supply tube | 50 | 1600-12 | 1650-12 |
| Adult cannula with 13' supply tube | 50 | 1600-13 | 1650-13 |
| Adult cannula with 14' supply tube | 50 | 1600-14 | 1650-14 |
| Adult cannula with 15' supply tube | 50 | 1600-15 | 1650-15 |
| Adult cannula with 16' supply tube | 50 | 1600-16 | 1650-16 |
| Adult cannula with 20' supply tube | 25 | 1600-20 | 1650-20 |
| Adult cannula with 21' supply tube | 25 | 1600-21 | 1650-21 |
| Adult cannula with 24' supply tube | 25 | 1600-24 | 1650-24 |
| Adult cannula with 25' supply tube | 25 | 1600-25 | 1650-25 |
| Adult cannula with 30' supply tube | 20 | 1600-30 | 1650-30 |
| Adult cannula with 35' supply tube | 20 | 1600-35 | 1650-35 |
| Adult cannula with 40' supply tube | 20 | 1600-40 | 1650-40 |
| Adult cannula with 50' supply tube | 20 | 1600-50 | 1650-50 |
| Adult cannula with 75' supply tube | 15 | 1600-75 | 1650-75 |
| Adult cannula with 100' supply tube | 10 | 1600-100 | 1650-100 |

Special lengths available on request. Minimum order is one (1) case.

Conventional Style Cannulas

| | Units | Safety Channel | Smooth Bore |
|---|-------|-------------------|-------------|
| Conventional cannula with 7' supply tube | 50 | 1053 ^x | 1056 |
| Conventional cannula without tube | 50 | 1057 ^x | |
| Conventional cannula with 16' supply tube | 50 | 1063 | 1068 |

* This cannula is available with closed-cell, over-the-ear, foam cushions, by adding TLC to the part number. (see catalog page SLC-00)



100 W. Sycamore Road - Arvin, CA 93203
 Phone: 661-854-3166
 1-800-421-0024 (CA only)
 1-800-235-4203 (U.S. & Canada)
 Fax: 661-854-3850
 Toll Free Fax: 1-800-628-4690 (U.S. & Canada)
 International Fax: 661-854-3850
 E-mail: salterlabs@us.salterlabs.com
 www.salterlabs.com



0482

ISO 9001/EN 46001 & ISO 13485