

Ventlab Corporation

AirFlow Resuscitation Bags with Manometer



- Integrated manometer promotes proper pressure delivery and thus aids in reducing gastric insufflations. This reduces vomiting, subsequent aspiration and other negative outcomes such as aspiration pneumonia.
- Vivid color coded circular manometer dial allows for easy and accurate determination of delivered pressures resulting in greater patient safety and clinician confidence.
- Available with disposable Bacterial/Viral filter attached to the exhalation port to help ensure an extra margin of safety for clinicians and staff.
- Optional 40 cm/H₂O or 25 cm/H₂O pop off valves help eliminate dangerously high delivered pressures for different patient populations.
- Nine different mask sizes from 00 to 7 help ensure a proper fit and the best possible ventilation technique for positive patient outcomes.
- Attached CO₂ indicators available on all bag sizes. Pre attachment of the CO₂ indicator helps ensure immediate use of the indicator, proper placement of the ET tube and the best possible ventilation practice.
- PEEP valves available on all bag sizes have an adjustable feature which allows for a wide range of settings to meet the varying needs of all patient populations.
- Custom configurations available on all sizes to help meet the specific needs of clinicians to aid in delivering the best possible patient care.

For more information or a product demonstration please contact your local Sales Representative (www.tri-anim.com).

Ordering Information

87-AF1140MB

Large Adult with Bag Reservoir, Manometer, Mask and Filter

87-AF1140MC

Large Adult with Corrugated Reservoir, Manometer, Mask and Filter

87-AF5140MB

Small Adult with Bag Reservoir, Manometer, Mask and Filter

87-AF5140MC

Small Adult with Corrugated Reservoir, Manometer, Mask and Filter

87-AF2140MB

Pediatric with Bag Reservoir, Manometer, Mask and Filter

87-AF2140MC

Pediatric with Corrugated Reservoir, Manometer, Mask and Filter

87-AF3140MB

Infant with Bag Reservoir, Manometer, Mask and Filter

87-AF3140MC

Infant with Corrugated Reservoir, Manometer, Mask and Filter

87-AF4140MB

Small Child with Bag Reservoir, Manometer, Mask and Filter

87-AF4140MC

Small Child with Corrugated Reservoir, Manometer, Mask and Filter

Custom configurations available.

Tri-anim is a specialty provider of medical services, supplying a wide variety of products paired with clinical expertise to support hospitals and other advanced patient-care facilities. With over 30 years of experience, Tri-anim understands the business of healthcare and offers innovative and value-oriented technologies, programs and modalities designed to improve patient care and maximize operational efficiencies.



800.874.2646 | www.tri-anim.com