

Ventlab Corporation

Hyperinflation Systems



- Integrated pressure manometer promotes proper pressure control which aids in the safest possible use of the product thus creating an extra margin of safety for the patient.
- Vivid color coded “circular” manometer dial allows for easy and accurate determination of pressures resulting in greater patient safety and clinician confidence.
- Design of the product puts both the manometer dial and the patient in the clinicians’ field of view at the same time. This allows for better assessment of the patient during the ventilation process.
- Patient port swivel allows for movement as necessary so the clinician can be in the best possible position to deliver optimum care.
- Thumb wheel dial is easy to use and activate. This results in the ability to fine tune pressures within the system.
- Multiple size masks are available to help ensure proper fit which promotes better patient outcomes.

Ordering Information

Use the table below to customize your product. Replace the first red number with the desired size face mask and the second red number for bag size as indicated in the example. (Example: 87-HS4**2**51 = Infant Mask and ½ liter Bag)

At the end of the product number, add **-M0** for the Small Neonate Mask or **-M00** for the Extra Small Neonate Mask. (Example: HS4**0**51**-M0** = Small Neonate Mask and 1/2 liter Bag)

87-HS4**XX**1

Series Hyperinflation Bag with Manometer

87-HS1**XX**0

Series Hyperinflation Bag w/o Manometer

Mask Sizes	Bag Sizes
0 = No Mask	5 = ½ liter
-M00 = X-Small Neonate	4 = ¼ liter
-M0 = Small Neonate	3 = 3 liter
1 = Small Infant	2 = 2 liter
2 = Infant	1 = 1 liter
3 = Small Child	
4 = Child	
5 = Adult (regular)	
6 = Adult (large)	
7 = Adult (reg. wide)	

Tri-anim is a specialty provider of medical services, supplying a wide variety of products paired with clinical expertise to support hospitals and other advanced patient-care facilities. With over 30 years of experience, Tri-anim understands the business of healthcare and offers innovative and value-oriented technologies, programs and modalities designed to improve patient care and maximize operational efficiencies.

For more information or a product demonstration please contact your local Sales Representative (www.tri-anim.com).

 **Tri-anim**[®]
800.874.2646 | www.tri-anim.com