



## LarySeal™ Pro


the **new generation** mask for safe and effective airway management

**LarySeal™ Pro** aids rapid and secure airway management, improves patient safety with gastric access to reduce risk of pulmonary aspiration, and speeds ET tube placement with a unique guidance system.




### Reduced Risk of Aspiration

Integrated suction catheter port for efficient removal of fluids and gastric contents.



### Directed Intubation

Accepting a range of ET tube diameters, the guide system directs the tip through the vocal cords.



### Flexible Tubing & Reinforced Bite Guard

Minimizing risk of pressure to soft tissues and protects from occlusion if patient bites down while also reducing likelihood of dental damage.

### Versatility

Standardizing on **LarySeal™ Pro** for both routine and advanced use means improved quality and safety, less training, and reduced costs with less stock holding. **LarySeal™ Pro** is the only mask you'll need.

## LarySeal™ Pro

**LarySeal™ Pro** features a **symmetrical cuff** and curvature that better suits the oropharyngeal anatomy for **quick** and **easy first time insertion** success.

### Ideal for:

- Routine and advanced airway management for patients with a higher risk of a difficult airway
- Being prepared for unexpected difficulties in intubation and airway control during anesthesia
- Airway rescue after failed intubation
- Outside OR airway management by experts and novices, especially during cardiac arrest
- Stomach decompression ability



Sterile ✓

Single Use ✓

DEHP Free ✓

Latex Free ✓



### Epiglottic Flap

Integral fenestrated flap protects from blockage with minimal increase in flow resistance, and elevates the epiglottis for easy, trauma-free ET tube or bronchoscope insertion.



### Clear Markings

Clear indicators with all the information a user needs including:

- Positioning and Depth
- Size
- Patient Weight
- Maximum ET Tube Diameter
- Maximum Suction Catheter Size
- Color-Coded Pilot Balloon

		PATIENT WEIGHT (KG)	SIZE	MAX ET TUBE ID (MM)	MAX SUCTION CATHETER (FR)	PACK SIZE
038-94-830U	LarySeal Pro Size 3	30kg - 50kg	3	7.5	16	10
038-94-840U	LarySeal Pro Size 4	50kg - 70kg	4	8	18	10
038-94-850U	LarySeal Pro Size 5	70kg - 100kg	5	8	18	10