



The Monaghan Advantage

- ❖ Global leaders in development and manufacture of aerosol drug delivery devices and respiratory management
- ❖ State-of-the-art aerosol research laboratory
- ❖ Patient focused features to promote ongoing patient compliance
- ❖ Efficacy supported by hundreds of peer reviewed articles and publications
- ❖ Focused on developing cost-efficient, outcomes-based solutions

ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
60050T	FILTERED MOUTHPIECE KIT Filter kits include a mouthpiece, one-way valved tee and a filter assembly. They are individually packaged with instructions.	cs/50

Contact your Monaghan Sales Territory Manager for more information or a demonstration.



REGISTERED
ISO 13485
QUALITY SYSTEM

Words or phrases accompanied by "™" and "®" are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation. © 2020 Monaghan Medical Corporation.

† trademarks and registered trademarks of the respective companies.

P/N 76364-00, 11/2020

Monaghan Medical
Corporation
Plattsburgh, NY 12901

518.561.7330 **phone**
800.833.9653 **toll free**
518.561.5088 **fax**

monaghanmed.com



FILTERED MOUTHPIECE KIT

- Filters exhaled emissions, viruses† and bacteria†.
- Reduces exposure of potentially harmful emissions to others in the room.
- Direct inhalation pathway ensures no loss of medication for the patient.
- No added resistance on exhalation.
- Tamper resistant valve.
- Standard 22 mm fittings to avoid unnecessary adapters.
- Single patient, disposable device.

FILTRATION PERFORMANCE

Viral Filtration Efficiency† >99.997%

Bacterial Filtration Efficiency† >99.9997%

†The bacterial and viral filtration efficiency performance of the filtration media has been tested by the material's manufacturer using a method adapted from ASTM F2101 in compliance with US FDA good manufacturing practice regulations 21 CFR Parts 210, 211 and 820.