




Curaplex® CuraBronch Elite™


Single-Use Bronchoscopes

Empowering clinical decisions
with clarity and control




Why Choose Curaplex® CuraBronch Elite™

CuraBronch Elite™ delivers high-definition imaging, precise control, and versatility to support caregivers in dynamic clinical environments.




Capture Photo and Video on Handle

Still images and video at the point of care, reducing the need to use monitor screen.




180° Scope Rotation Knob

Reduces the need to maneuver hand/arm position while holding the scope.



Swivel Suction Connector

Adjust the direction/position of suction tubing based on the user's unique setup.



Biopsy Valve

Standard offering with both slip tip and luer lock syringe connections.

High-Definition Imaging

Wide field of view (120°), 720 x 720 camera chip resolution and an HD 1080p monitor screen.





Built to Support a Range of Bronchoscopy Procedures

Large Working Channel

Equal to or larger working channel across all three sizes compared to market-leading bronchoscopes of the same size.¹

Slim Insertion Tube Profile


Equal to or smaller insertion tube and smaller distal tip outer diameter than market-leading bronchoscopes of the same size.¹

±180° Tip Flexion

Navigate complex anatomy precisely across all three sizes.

ET Tube Retainer


Built-in ET tube retainer for use during intubation, along with an Aintree Intubation Catheter retainer in the slim size bronchoscope.



1 - Competitive comparison based on publicly available technical data sheets for Ambu® aScope™ 4 Broncho and Verathon BFlex™ 2 single-use bronchoscopes. Like-for-like sizes compared (Slim/Regular/Large).

Full-HD Display for Real-Time Decision Support

- High definition 1080p resolution on a 14" screen with brightness, contrast, sharpness, and color controls
- Up to 5.5 hours of battery life on a single charge
- Connect to external monitor with the HDMI port, export data via USB (DICOM format supported)
- Can be mounted to workstation, universal VESA pole mount, or propped on tabletop with kickstand



Ordering Information

Select the right scope for the right procedure



Single-Use Bronchoscopes

Item #	Product	Qty
301-ECB-100S	CuraBronch Elite™ Slim (3.7/1.2)	5
301-ECB-100R	CuraBronch Elite™ Regular (4.9/2.2)	5
301-ECB-100L	CuraBronch Elite™ Slim (5.7/3.0)	5

Monitor

Item #	Product	Qty
301-EDS-500	Curaplex® HD Video Monitor	1

Workstation

Item #	Product	Qty
301-ECA-200	Curaplex® Elite Workstation	1

More Information

tri-anim.com | 800.874.2646

or contact your Tri-anim® account manager

Portable and Compact Workstation that Completes the Curabronch Elite™ Package

- Compatible with Curaplex® HD Video Monitor
- Adjustable arm for direct positioning in your line of sight
- Large prep tray for easy access to airway tools
- Scope holder plus storage hook for unopened scopes
- Storage bin for additional organization



For prescription use only. Caution: Federal law (U.S.) restricts this device to sale by or on the order of a physician. The Curaplex® CuraBronch Elite™ is intended for use in a hospital environment or medical office environment. It is designed for use in adults.

The Curaplex® CuraBronch Elite™ is a sterile, single-use, and flexible device intended to be operated with its compatible display system Curaplex® HD Video Monitor (301-EDS-500). The device provides endoscopic visualization of airways and tracheobronchial tree, in conjunction with displaying system to perform various endoscopic procedures.