



VibraPEP®



The VibraPEP® by Curaplex® is a single-patient multiple use device that provides oscillatory positive expiratory pressure (OPEP) therapy for the detachment and removal of pulmonary secretions. Through five adjustable therapy settings, clinicians can tailor expiratory resistance and oscillation to the individual patient needs. The VibraPEP® uses the patient's full expired air volume to produce pressure and oscillatory vibrations. The success of OPEP therapy depends on these vibrations, especially for patients with a low expiratory volume.

- CombinedPEP® promotes collateral ventilation through positive expiratory pressure and DynamicPEP® generates asynchronous flow and pressure oscillations that assist in the detachment and removal of bronchial secretions.
- Oscillations begin at flows as low as 10 lpm.
- Common ISO fittings permit versatility of patient interfaces.
- Unique "valve tube" design is not gravity dependent and can be used in any position.

Item #	Description	Qty
301-44F50-100	VibraPEP® with Valved Tee Adapter	10/bx
301-44F10-120	VibraPEP® without Valved Tee Adapter	10/bx
301-12218	Valved Tee Adapter	10/bx

VibraPEP®

