

Telesair® Bonhawa High Flow Oxygen Therapy System

Bon (French: “good”) +
hawa (Hawaiian: “home”)

The next generation in High Flow
Oxygen Therapy, designed with clinicians
in mind and now exclusively* available
through Tri-anim Health Services.



*not exclusive in government accounts

Bonhawa is a registered trademark of Telesair, Inc. © 2025
Tri-anim Health Services. All rights reserved.

High Flow Oxygen Therapy You Can Confirm in Real Time

Bonhawa represents the next generation of high flow therapy, engineered around clinician needs with faster disinfection processing time, flows up to 70 LPM, and Confirmed Therapy™ – real-time flow monitoring that automatically compensates for pressure fluctuations so clinicians can see exactly what the patient receives. Drive better outcomes with therapeutic humidification across all flows and quick therapy set-up – no delays, just optimal patient care. Bonhawa enables more time for patient care and less time on device management.



PERFORMANCE

Advanced respiratory support when it matters most

Extended High Flow Oxygen Therapy (HFOT) range up to 70 LPM

Confirmed Therapy - real-time flow monitoring that automatically compensates for pressure fluctuations, confirming prescribed therapy is delivered

Maintain therapeutic humidification without disruption¹




EFFICIENCY

Streamlined workflow, reduced cost of ownership

5 minutes or less disinfection time provides labor and cost savings

Intuitive forward facing interface for fast, confident setup

Fast setup and minimal maintenance with exterior gas pathway



VERSATILITY

Adaptable solutions for your clinical scenario

Full-capability transport - don't sacrifice flow, heat, and humidity for up to 90 minutes (when used with UPS)

Ability to maintain heat and humidity parameters on adult and pediatric patients during transport

Universal cannulas and med delivery compatibility

¹ when used with UPS Battery both heat and humidity can be maintained during transport up to 90 minutes

BUILD YOUR BONHAWA HIGH FLOW OXYGEN THERAPY SOLUTION

*Customizable to your needs—
connect with your account manager to
create the right configuration for you.*

BONHAWA DEVICE

- Bonhawa Respiratory Humidifier (P/N: 7147-BHW01)
- 2-year warranty
- Serviceable and upgradeable

CONSUMABLES

- Heated Breathing Circuit Pack (P/N: 7147-HBK01-10)
- Nasal Cannulas
 - Large (P/N: 7147-ENC01-10)
 - Medium (P/N: 7147-ENC02-10)
 - Small (P/N: 7147-ENC03-10)
- Air Filter (P/N: 7147-AFK01-10)

ACCESSORIES

- Zopec Transport Battery (P/N: 7147-UPS01)
- Medical Trolley with Bonhawa Tray and Basket (P/N: 7147-CTP02)
- Ventilator Support Arm (P/N: 7147-VSA01)
- Oxygen Tank Holder (P/N: 7147-OCH01)
- Precision Manifolds - Pediatric/Adult/Yblock etc (715-509080)



ORDERING INFORMATION

Setup requirements include the Telesair Bonhawa HFOT Humidifier, Breathing Circuit Disposable Kit, 22mm Nasal Cannula of choice, and 70LPM Flowmeter Manifold setup.

Item #	Description	Qty
7147-BHW01	Bonhawa Respiratory Humidifier	1/ea
7147-HBK01-10	Heater Breathing Circuit Set with water chamber adapter, a water chamber, and a heated breathing circuit	10/pk
715-509080	Telesair Flowmeter and Oxygen set up includes: Quad block with O2 connection, male one-way valves, 1 70LPM flowmeter, 2 3-ft O2 hoses, and 1-9ft O2 hose, DISS connections	1/bx

Transport setup requirements include the Medical Trolley with basket included and the Oxygen Tank Holder.

Item #	Description	Qty
7147-UPS01	Zopec UPS Transport Battery	1/ea
7147-OCH01	Bonhawa Oxygen Tank Holder, holds 2 transport tanks	1/ea
7147-CTP02	Medical Trolley (3 pieces, Aluminum/SS) including base with 5 castors/2 with brake, Trolley Pole and Basket	1/ea

The high flow nasal cannulas below pair seamlessly with Bonhawa for high flows up to 70LPM.

Standard 22mm connection. Choose from the below cannula sizing to meet your needs:

Item #	Description	Size	Application	Flow Range	Length	Qty
7147-ENC03-10	Telesair High Flow Nasal Cannula	Small	Pediatric or Adult Low Flow	2-40 LPM	32 cm	10/pk
7147-ENC02-10	Telesair High Flow Nasal Cannula	Medium	Adult Standard	10-80 LPM*	32 cm	10/pk
7147-ENC01-10	Telesair High Flow Nasal Cannula	Large	Adult High-Flow	20-80 LPM*	35 cm	10/pk

Additional Accessories:

Ventilator Support Arm helps provide medication delivery with proper setup on the Bonhawa device. The Air Filter is recommended to change every 30 days on the Bonhawa device.

Item #	Description	Qty
7147-VSA01	Ventilator Support Arm	1/ea
7147-AFK01-10	Bonhawa Air Filters for intake air, blended synthetic fiber	30/bx
715-509080	Precision Manifolds - Pediatric/Adult/Yblock etc	1/ea

TELSAIR

Tri-anim