

Telesair® Bonhawa **High Flow Oxygen** **Therapy System**

Bon (French: “good”) +
hawa (Hawaiian: “home”)

The next generation in High Flow Oxygen Therapy, designed with clinicians in mind and now exclusively* available through Tri-anim Health Services.



*not exclusive in government accounts

Bonhawa is a registered trademark of Telesair, Inc. © 2025
Tri-anim Health Services. All rights reserved.

Telesair® Bonhawa

Bonhawa represents the next generation of high flow therapy, engineered around clinician needs with faster disinfection processing time, flows up to 70LPM, and real-time flow monitoring that also monitors temperature and oxygen through our large front-facing touchscreen. Drive better outcomes with therapeutic humidification across all flows and quick therapy set-up - no delays, just optimal patient care. Bonhawa enables more time for patient care and less time on device management.



PERFORMANCE

***Advanced respiratory support
when it matters most***

Extended High Flow
Oxygen Therapy (HFOT)
range up to 70 LPM

Real-time flow monitoring
confirms therapy

Maintain therapeutic
humidification
without disruption¹



EFFICIENCY

***Streamlined workflow,
reduced cost of ownership***

5 minutes or less
disinfection time provides
labor and cost savings

Intuitive interface,
forward facing screen,
quickly set up patients

Fast setup and
minimal maintenance



VERSATILITY

***Adaptable solutions for your
clinical scenario***

Continuous therapy
(up to 90 min) on-the-go¹

Ability to maintain heat
and humidity parameters
on adult and pediatric
patients during transport

Universal cannulas and
med delivery compatibility

¹ when used with UPS Battery both heat and humidity can be maintained during transport up to 90 minutes

BUILD YOUR BONHAWA HIGH FLOW OXYGEN THERAPY SOLUTION

BONHAWA DEVICE

- Bonhawa Respiratory Humidifier (P/N: 7147-BHW01)
- 2-year warranty
- Serviceable and upgradeable



CONSUMABLES

- Heated Breathing Circuit Pack (P/N: 7147-HBK01-10)
- Nasal Cannulas
 - Large (P/N: 7147-ENC01-10)
 - Medium (P/N: 7147-ENC02-10)
 - Small (P/N: 7147-ENC03-10)
- Air Filter (P/N: 7147-AFK01-10)



ACCESSORIES

- Zopec Transport Battery (P/N: 7147-UPS01)
- Medical Trolley with Bonhawa Tray and Basket (P/N: 7147-CTP02)
- Ventilator Support Arm (P/N: 7147-VSA01)
- Oxygen Tank Holder (P/N: 7147-OCH01)



ORDERING INFORMATION

Setup requirements include the Telesair Bonhawa HFOT Humidifier, Breathing Circuit Disposable Kit, 22mm Nasal Cannula of choice, and 70LPM Flowmeter Manifold setup.

Item #	Description	Qty
7147-BHW01	Bonhawa Respiratory Humidifier	1/ea
7147-HBK01-10	Heater Breathing Circuit Set with water chamber adapter, a water chamber, and a heated breathing circuit	10/pk
715-509080	Telesair Flowmeter and Oxygen set up includes: Quad block with O2 connection, male one-way valves, 1 70LPM flowmeter, 2 3-ft O2 hoses, and 1-9ft O2 hose, DISS connections	1/bx

Transport setup requirements include the Medical Trolley with basket included and the Oxygen Tank Holder.

Item #	Description	Qty
7147-UPS01	Zopec UPS Transport Battery	1/ea
7147-UCH01	Bonhawa Oxygen Tank Holder, holds 2 transport tanks	1/ea
7147-CTP02	Medical Trolley (3 pieces, Aluminum/SS) including base with 5 castors/2 with brake, Trolley Pole and Basket	1/ea

The high flow nasal cannulas below pair seamlessly with Bonhawa for high flows up to 70LPM.

Standard 22mm connection. Choose from the below cannula sizing to meet your needs:

Item #	Description	Size	Application	Flow Range	Length	Qty
7147-ENC03-10	Telesair High Flow Nasal Cannula	Small	Pediatric or Adult Low Flow	2-40 LPM	32 cm	10/pk
7147-ENC02-10	Telesair High Flow Nasal Cannula	Medium	Adult Standard	10-80 LPM*	32 cm	10/pk
7147-ENC01-10	Telesair High Flow Nasal Cannula	Large	Adult High-Flow	20-80 LPM*	35 cm	10/pk

Additional Accessories:

Ventilator Support Arm helps provide medication delivery with proper setup on the Bonhawa device. The Air Filter is recommended to change every 30 days on the Bonhawa device.

Item #	Description	Qty
7147-VSA01	Ventilator Support Arm	1/ea
7147-AFK01-10	Bonhawa Air Filters for intake air, blended synthetic fiber	30/bx

TELSAIR

Tri-anim